



Ball Abutment

Specifications

- Use to create overdenture for edentulous patient
- Path recovery to the maximum of 15°
- Tightening torque: 30 Ncm
- Use ball abutment driver

Mini Hex Connections						
	GH	1.0	2.0	3.0	4.0	5.0
	REFERENCE CODE	ASC1M	ASC2M	ASC3M	ASC4M	ASC5M

Regular Hex Connections						
	GH	1.0	2.0	3.0	4.0	5.0
	REFERENCE CODE	ASC1	ASC2	ASC3	ASC4	ASC5