


# Temporary & Analog


## Specifications

- For temporary crown on single or multi implant
- Tightening Torque: 30 Ncm
- Use 1.2 hex screw driver

## Temporary Abutment


Mini Hex Connections		
	DIAMETER	Hex      Non-Hex
	GH	1      1
	REFERENCE CODE	TA400HM      TA400NM


Abutment Screw included

Regular Hex Connections		
	DIAMETER	Hex      Non-Hex
	GH	1      1
	REFERENCE CODE	TA450H      TA450N

Abutment Screw included

## Analog

Mini Hex Connections		
	DIAMETER	3.0 , 3.5
	REFERENCE CODE	RL300M

Regular Hex Connections			
	DIAMETER	3.7	4.0 ,4.5 ,5.0 ,6.0
	REFERENCE CODE	RL350	RL450