


Impression Coping

Closed Tray Impression Coping (Pop In)

- For impression on single implant
- For impression on multiple implants when all are parallel
- Tightening torque : 5-10 Ncm
- Use 1.2 hex screw driver

Mini Hex Connections					
	DIAMETER (ϕ)	4.0	4.0	4.5	4.5
	HEIGHT	SHORT	LONG	SHORT	LONG
	REFERENCE CODE	TI400M	TI400LM	TI450M	TI450LM

Regular Hex Connections							
	DIAMETER (ϕ)	4.5	4.5	5.5	5.5	6.5	6.5
	HEIGHT	SHORT	LONG	SHORT	LONG	SHORT	LONG
REFERENCE CODE	TI450	TI450L	TI550	TI550L	TI650	TI650L	